# NADAL® CRP Test (test cassette)

REF 311801N-10/311801N-20

|           |                            |    |           |                    |    |
|-----------|----------------------------|----|-----------|--------------------|----|
| <b>DE</b> | Gebrauchsanweisung         | 2  | <b>CZ</b> | Návod k použití    | 27 |
| <b>EN</b> | Instructions for use       | 7  | <b>FI</b> | Käyttöohje         | 31 |
| <b>FR</b> | Instructions d'utilisation | 11 | <b>NL</b> | Gebruiksaanwijzing | 35 |
| <b>ES</b> | Instrucciones de uso       | 15 |           | Symbols            | 43 |
| <b>IT</b> | Istruzioni per l'uso       | 19 |           | Our Teams          | 44 |
| <b>PL</b> | Sposób użycia              | 23 |           |                    |    |

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### 16. Referenzen

1. Morley JJ, Kushner (1982) Serum C-reactive protein levels in disease. In: Kushner I, Volanakis JE, Gewurz H, eds. C-reactive protein and the plasma protein response to tissue injury. Ann. NY Acad. Sci. 389: 406-417
2. Van Lente F (1982) The Diagnostic Utility of C-Reactive Protein. Hum Path 13(12): 1061-3.
3. Peltola HO (1982) C-reactive protein for rapid monitoring of infections of the central nervous system. Lancet:980-983.
4. Shaw AC (1991) Serum C-Reactive Protein and Neopterin Concentrations in Patients with Viral or Bacterial Infection. J Clin Pathol 44(7): 596-9.
5. Dowton SR and Colten HR (1988) Acute Phase Reactants in Inflammation and Infection. Semin Hematol 25(2): 84-90.
6. Gambino R (1994) C-Reactive Protein (CRP) - How much Proof do we need? Lab Rep 16(11): 83-5
7. Thomas L (2012) Labor und Diagnose – Indikationen und Bewertung von Laborbefunden für die medizinische Diagnostik. TH-Books Verlagsgesellschaft mbH.
8. Auflage. Band 2. S.1279.

Rev.0, 2020-01-28 OM

### 1. Intended Use

The NADAL® CRP Test is a lateral flow chromatographic immunoassay for the semiquantitative detection of C-reactive protein (CRP) in human whole blood, serum or plasma specimens. The test is intended for use as an aid in the diagnosis of bacterial infectious diseases and inflammatory processes (see section 12 'Limitations'). Due to a large number of possible symptoms, the test is not restricted to a defined target patient group and can generally be used to differentiate between bacterial and viral infections or to assess the presence or severity of inflammatory processes. The test procedure is not automated and requires no special training or qualification. The NADAL® CRP Test is designed for professional use only.

### 2. Introduction and Clinical Significance

C-reactive protein (CRP) is detected in patient specimens in conjunction with acute infections, necrosis and various inflammatory disorders.<sup>1,2</sup> CRP values do not increase to the same extent in viral infections as in bacterial ones, which is why this parameter can contribute to their differentiation.<sup>3,4</sup> CRP is a protein which is predominantly formed by the liver during an acute inflammatory process. A positive test result indicates the presence, not the cause, of an acute inflammatory reaction. CRP synthesis is stimulated by antigen-immune complexes, bacteria, fungi or trauma. Monitoring of CRP levels in patients can provide information on the effectiveness of treatment and facilitate the assessment of patient recovery.<sup>5</sup>

The inter-individual variation of CRP values is relatively high. Age, overweight, diabetes mellitus type 2, smoking, intake of gestagen (progestin)/estrogen medication or atherosclerotic/cardiovascular diseases, among others, may contribute to elevated values. In general, values of  $>10 \mu\text{g/mL}$  can be regarded as elevated for the majority of patients.

Disease severity can be classified on the basis of the CRP concentration into:

- Mild inflammation ( $>10\text{--}40 \text{ mg/L}$ ). Possible causes: local abscess, mild operative or accident trauma, heart attack, deep vein thrombosis, inactive rheumatic diseases, metastasised malignant tumour and isolated viral infections.
- Moderate inflammation ( $>40\text{--}100 \text{ mg/L}$ ). Possible causes: severe inflammatory processes like purulent cystitis, bronchitis, dental suppuration, urinary tract infections and genital infections.
- High-grade inflammation ( $>100 \text{ mg/L}$ ). Possible causes: acute generalised bacterial or fungal infections (sepsis) and severe tissue injury following polytrauma or major surgical procedures.<sup>7</sup>

### 3. Test Principle

The NADAL® CRP Test is a lateral flow chromatographic immunoassay for the semiquantitative detection of CRP in human whole blood, serum or plasma specimens.

After the specimen is added to the specimen well (S) of the test cassette, it reacts with anti-CRP antibodies (active ingredient of the detection reaction: approx.  $0.1 \mu\text{g/test}$ ) which are conjugated to coloured particles and precoated onto the conjugate pad of the internal test strip. The mixture

then migrates along the membrane chromatographically by capillary action and reacts with further anti-CRP antibodies (active ingredient of the detection reaction: approx.  $0.5 \mu\text{g/test}$ ) immobilised in the test line regions: 'T1', 'T2' and 'T3' of the membrane. The presence of a coloured line in the test line region(s) 'T1' and/or 'T2', and/or 'T3' indicates a positive result. The absence of a coloured line in the test line region(s) 'T1' and/or 'T2', and/or 'T3' indicates a negative result. The number of lines depends on the CRP concentration in the specimen. The higher the CRP concentration in the specimen, the more coloured lines become visible.

The formation of a coloured line in the control line region 'C' serves as a procedural control, indicating that the proper volume of specimen has been added and membrane wicking has occurred.

### 4. Reagents and Materials Supplied

- 10/20 NADAL® CRP test cassettes\*
  - 10/20 disposable pipettes
  - 10/20 end-to-end capillary tubes ( $10 \mu\text{L}$ )
  - 10/20 extraction tubes containing buffer\*\*
  - 1 package insert
- \*containing the preservative sodium azide  
 \*\*Phosphate buffer containing the following preservatives:  
 Kanamycin sulfate:  $0.25 \text{ g/L}$  and ProClin™ 300:  $<0.03\%$ .

No hazard labelling is required according to Regulation (EC) № 1272/2008 CLP. Concentrations are below exemption threshold.

### 5. Additional Materials Required

- Specimen collection containers (appropriate for specimen material to be tested)
- Centrifuge (for serum or plasma specimens only)
- Alcohol pads
- Lancets (for fingerstick whole blood specimens only)
- Timer

### 6. Storage & Stability

Test kits should be stored at  $2\text{--}30^\circ\text{C}$  until the indicated expiry date. Test cassettes are stable until the expiry date printed on the foil pouches. Test cassettes must remain in the sealed foil pouches until use. Do not freeze the test kit. Do not use tests beyond the expiry date indicated on the packaging. Care should be taken to protect test kit components from contamination. Do not use test kit components if there is evidence of microbial contamination or precipitation. Biological contamination of dispensing equipment, containers or reagents can lead to inaccurate results.

### 7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Carefully read through the test procedure prior to testing.
- Do not use the test beyond the expiration date indicated on the packaging.
- Do not use test kit components if the primary packaging is damaged.
- Tests are for single use only.
- Do not add specimens to the reaction area (result area).
- In order to avoid contamination, do not touch the reaction area (result area).

- Avoid cross-contamination of specimens by using a new extraction tube for each specimen obtained.
- Do not substitute or mix components from different test kits.
- Do not eat, drink or smoke in the area where specimens and test kits are handled.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being assayed.
- Handle all specimens as if they contain infectious agents. Observe established precautions for microbiological risks throughout all procedures and standard guidelines for the appropriate disposal of specimens.
- The test kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled in accordance with usual safety precautions (e.g., do not ingest or inhale).
- Temperature can adversely affect test results.
- Used testing materials should be disposed of according to local regulations.

## 8. Specimen Collection and Preparation

The NADAL® CRP Test can be performed using whole blood (from venipuncture or fingerstick), serum or plasma.

### To collect fingerstick whole blood specimens:

- Wash the patient's hand with soap and warm water or clean it with an alcohol pad. Allow it to dry.
- Massage the hand, without touching the puncture site, by rubbing along the hand towards the fingertip of the middle or ring finger.
- Puncture the skin with a sterile lancet. Wipe away the first drop of blood.
- Gently rub the hand from the wrist to the palm, and then to the finger to form a rounded drop of blood over the puncture site.

Fingerstick whole blood should be tested immediately.

### Venipuncture whole blood specimens

Containers containing anticoagulants, such as K<sub>2</sub>EDTA, sodium citrate, sodium heparin or lithium heparin should be used for the preparation of venous whole blood or plasma specimens.

Testing should be performed immediately after specimen collection. Do not leave specimens at room temperature for prolonged periods of time.

If the test is to be run within 2 days of specimen collection, whole blood collected by venipuncture should be stored at 2-8°C.

Do not freeze whole blood specimens.

### Serum and plasma specimens

Separate serum or plasma from blood as soon as possible to avoid haemolysis. Use only clear, non-haemolysed specimens.

Testing should be performed immediately after specimen collection. Do not leave specimens at room temperature for prolonged periods of time. Serum and plasma specimens can

be stored at 2-8°C for up to 3 days. For long-term storage, specimens should be kept at -20°C.

Bring specimens to room temperature prior to testing. Frozen specimens should be completely thawed and mixed well prior to testing. Specimens should not be frozen and thawed repeatedly.


If specimens are to be shipped, they should be packed in compliance with all applicable regulations for the transportation of etiologic agents.

Icteric, lipemic, haemolysed, heat-treated and contaminated specimens may lead to inaccurate test results.

## 9. Test Procedure

**Bring tests, specimens, buffer and/or controls to room temperature (15-30°C) prior to testing.**

1. Unscrew the extraction tube's green cap.


2. **a) For serum or plasma specimens:**

Holding the pipette vertically, draw the specimen up to the first widening (approximately 5 µL) and add it to the extraction tube containing buffer.

**b) For venipuncture whole blood specimens:**


Holding the pipette vertically, add 1 drop (approximately 10 µL) of a whole blood specimen to the extraction tube containing buffer.

**c) For fingerstick whole blood specimens:**

Apply the end of the end-to-end capillary tube to the blood until it is filled completely (approximately 10 µL). **Avoid air bubbles forming.**

Place the blood-filled capillary into the extraction tube containing buffer.

3. Close the tube and invert it several times so that the specimen and the buffer get mixed well.


4. Remove the test cassette from the foil pouch and use it as soon as possible. The best results will be obtained if the test is performed immediately after opening the foil pouch. Label the test cassette with the patient or control identification.

5. Place the test cassette on a clean and level surface.

6. Unscrew the extraction tube's white cap.

7. Holding the extraction tube vertically, add 3 drops of the extracted specimen to the specimen well (S) of the test cassette.

**Avoid trapping air bubbles in the specimen well (S) and do not add any solution to the result area.**

8. Start the timer.
9. Wait for the coloured line(s) to appear. Read the test result after 5 minutes. Do not interpret the result after more than 8 minutes.


## 10. Result Interpretation

### Positive


#### CRP concentration: $\geq 10$ mg/L

A coloured line develops in the control line region 'C' and another coloured line develops in the test line region 'T1'. The result indicates that a CRP concentration of  $\geq 10$  mg/L has been detected in the specimen.


#### CRP concentration: $\geq 40$ mg/L

In addition to the control line 'C', a coloured line develops in the test line region 'T1' and another in the test line region 'T2'. The result indicates that a CRP concentration of  $\geq 40$  mg/L has been detected in the specimen.


#### CRP concentration: $\geq 80$ mg/L

In addition to the control line 'C', a coloured line develops in each of the test line regions 'T1', 'T2' and 'T3'. The result indicates that a CRP concentration of  $\geq 80$  mg/L has been detected in the specimen.

### Negative

A coloured line develops in the control line region 'C'. No line develops in the test line region 'T1', 'T2' or 'T3'. The result indicates that a CRP concentration of  $< 10$  mg/L has been detected in the specimen.

### Invalid

The control line 'C' fails to appear. Results from any test which has not produced a control line at the specified reading time must be discarded. Please review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your distributor.

Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for the control line failure.

## 11. Quality Control

An internal procedural control is included in the test cassette:

A coloured line appearing in the control line region 'C' is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

*Good laboratory practice (GLP)* recommends the use of external control materials to ensure proper test kit performance.

## 12. Limitations

- The NADAL® CRP Test is for professional *in-vitro* diagnostic use only. It should be used for the semiquantitative detection of CRP in human whole blood, serum or plasma specimens only.
- The NADAL® CRP Test only semiquantitatively determines a CRP level in human whole blood, serum and plasma specimens, and should not be used as the sole criterion for a diagnosis or for evaluating the risks of inflammatory disorders.
- As with all diagnostic tests, a definitive clinical diagnosis should only be made by a physician following the evaluation of all clinical and laboratory findings.
- Even in the case of a positive result, a comprehensive analysis of clinical symptoms and the use of other detection methods should always be considered.
- High titres of rheumatoid factors (RF) in positive clinical specimens may interfere with test results.
- There is a slight possibility that whole blood specimens with a very high viscosity or whole blood specimens which were stored for more than 2 days might cause migration problems. Using a new test cassette, repeat the test with a serum or plasma specimen from the same patient.
- Elevated levels of CRP in oral contraceptive (OC) users should be considered with caution as the American Physiological Society recommends further studies on the impact of OC use on CRP and other inflammatory parameters.
- High concentrations of CRP may produce a high dose hook effect, resulting in false negative results. A high dose hook effect has not been observed with the NADAL® CRP Test up to a CRP concentration of 3120 mg/L.
- As the test provides semiquantitative results only, the concentration gradient during monitoring can only be determined very roughly within the ranges listed above. If smaller concentration changes need to be detected, a quantitative test must be used.

## 13. Expected values

CRP is a non-specific marker for acute infections, necrotic conditions and a variety of inflammatory disorders. According to American Heart Association (AHA) guidelines, the expected CRP value should be  $< 10$  mg/L to rule out bacterial infections risks.

## 14. Performance Characteristics

### Clinical performance

#### Diagnostic sensitivity and specificity

The NADAL® CRP Test was evaluated using clinical specimens in comparison with a leading, commercially available CRP EIA.

| Method              |                | CRP EIA       |                |                |             |       |
|---------------------|----------------|---------------|----------------|----------------|-------------|-------|
|                     |                | Nega-<br>tive | Positive       |                |             | Total |
| NADAL® CRP Test     | Result         | 0-<10<br>mg/L | 10-<40<br>mg/L | 40-<80<br>mg/L | ≥80<br>mg/L |       |
|                     | 0-<10<br>mg/L  | 99            | 1              | 0              | 0           | 100   |
|                     | 10-<40<br>mg/L | 3             | 52             | 1              | 0           | 56    |
|                     | 40-<80<br>mg/L | 0             | 3              | 25             | 5           | 33    |
|                     | ≥80<br>mg/L    | 0             | 0              | 3              | 79          | 82    |
| Total               |                | 102           | 56             | 29             | 84          | 271   |
| % Overall agreement |                | 97.1%         | 92.9%          | 86.2%          | 94.0%       | 98.5% |
| 99.4%               |                |               |                |                |             |       |

Diagnostic sensitivity: 168/169=99.4% (CI\*: 96.7% - 100%)

Diagnostic specificity: 99/(99+3)=97.1% (CI\*: 93.9% - 98.1%)

Overall agreement: (99+168)/(99+3+169) = 98.5% (CI\*: 91.6%- 99.4%)

\*95% confidence interval

#### Analytical performance

##### Analytical sensitivity

Each of the specimens with CRP concentrations of 0 mg/L, 3 mg/L, 10 mg/L, 40 mg/L, 80 mg/L and 200 mg/L was tested in replicates of 5 using 3 lots of the NADAL® CRP Test. All specimens with a CRP concentration of ≥10 mg/L showed positive results.

##### Detection limit

The detection limit of the NADAL® CRP Test is 10 mg/L.

##### Interfering substances

Negative plasma and serum pools (confirmed by ELISA) as well as 10 mg/L CRP positive plasma and serum specimens (confirmed by ELISA) spiked with the following interfering substances showed no interference with the NADAL® CRP Test.

|                      |            |
|----------------------|------------|
| Acetaminophen        | 200 mg/L   |
| Acetylsalicylic acid | 200 mg/L   |
| Albumin              | 20000 mg/L |
| Ascorbic acid        | 200 mg/L   |
| Bilirubin            | 10000 mg/L |
| Caffeine             | 200 mg/L   |
| Creatine             | 2000 mg/L  |
| Gentisic acid        | 200 mg/L   |
| Haemoglobin          | 10000 mg/L |
| Methanol             | 10%        |
| Oxalic acid          | 600 mg/L   |

##### Cross-reactivity

HAV, HBsAg, HIV, HCV, HEV, syphilis, HAMA, RF, MONO, Rubella, TOXO, CMV, AFP and CEA positive specimens (confirmed by ELISA and a clinical diagnosis) and the same specimens but spiked with 20 mg/mL CRP at the dilution rate of 1:1 were tested using the NADAL® CRP Test.

No cross-reactivity with the specimens was observed when tested using the NADAL® CRP Test.

#### Measuring range

A high positive CRP specimen (3120 mg/L) was diluted with the same negative plasma to concentrations of 1:2, 1:16, 1:64, 1:128, 1:300, 1:600, 1:1000 and tested in triplicate using the NADAL® CRP Test. No high dose hook effect was observed when tested using the NADAL® CRP Test.

#### Precision

##### Repeatability

Repeatability was established by testing 10 replicates of 4 CRP specimens (0 mg/L, 10 mg/L, 40 mg/L and 80 mg/L).

##### Reproducibility

Reproducibility was established by testing 4 CRP specimens (0 mg/L, 10 mg/L, 40 mg/L and 80 mg/L) using 3 independent lots of the NADAL® CRP Test.

The NADAL® CRP Test demonstrated acceptable repeatability and reproducibility. The negative and positive values were correctly identified >99% of the time.

#### 15. Serious incident reporting

In the case of any serious incidents related to the performance of the NADAL® CRP Test, please inform nal von minden GmbH and the competent authority immediately. If still possible, do **not** dispose of the used test and the corresponding test kit components.

#### 16. References

- Morley JJ, Kushner (1982) Serum C-reactive protein levels in disease. In: Kushner I, Volanakis JE, Gewurz H, eds. C-reactive protein and the plasma protein response to tissue injury. Ann. NY Acad. Sci. 389: 406-417
- Van Lente F (1982) The Diagnostic Utility of C-Reactive Protein. Hum Path 13(12): 1061-3.
- Peltola HO (1982) C-reactive protein for rapid monitoring of infections of the central nervous system. Lancet:980-983.
- Shaw AC (1991) Serum C-Reactive Protein and Neopterin Concentrations in Patients with Viral or Bacterial Infection. J Clin Pathol 44(7): 596-9.
- Dowton SR and Colten HR (1988) Acute Phase Reactants in Inflammation and Infection. Semin Hematol 25(2): 84-90.
- Gambino R (1994) C-Reactive Protein (CRP) - How much Proof do we need? Lab Rep 16(11): 83-5
- Thomas L (2012) Labor und Diagnose – Indikationen und Bewertung von Laborbefunden für die medizinische Diagnostik. TH-Books Verlagsgesellschaft mbH. 8. Auflage. Band 2. S.1279.

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### 16. Bibliographie

- Morley JJ, Kushner (1982) Serum C-reactive protein levels in disease. In: Kushner I, Volanakis JE, Gewurz H, eds. C-reactive protein and the plasma protein response to tissue injury. Ann. NY Acad. Sci. 389: 406-417
- Van Lente F (1982) The Diagnostic Utility of C-Reactive Protein. Hum Path 13(12): 1061-3.
- Peltola HO (1982) C-reactive protein for rapid monitoring of infections of the central nervous system. Lancet:980-983.
- Shaw AC (1991) Serum C-Reactive Protein and Neopterin Concentrations in Patients with Viral or Bacterial Infection. J Clin Pathol 44(7): 596-9.
- Downton SR and Colten HR (1988) Acute Phase Reactants in Inflammation and Infection. Semin Hematol 25(2): 84-90.
- Gambino R (1994) C-Reactive Protein (CRP) - How much Proof do we need? Lab Rep 16(11): 83-5
- Thomas L (2012) Labor und Diagnose – Indikationen und Bewertung von Laborbefunden für die medizinische Diagnostik. TH-Books Verlagsgesellschaft mbH. 8. Auflage. Band 2. S.1279.

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### 16. Referencias

1. Morley JJ, Kushner (1982) Serum C-reactive protein levels in disease. In: Kushner I, Volanakis JE, Gewurz H, eds. C-reactive protein and the plasma protein response to tissue injury. Ann. NY Acad. Sci. 389: 406-417
2. Van Lente F (1982) The Diagnostic Utility of C-Reactive Protein. Hum Path 13(12): 1061-3.
3. Peltola HO (1982) C-reactive protein for rapid monitoring of infections of the central nervous system. Lancet:980-983.
4. Shaw AC (1991) Serum C-Reactive Protein and Neopterin Concentrations in Patients with Viral or Bacterial Infection. J Clin Pathol 44(7): 596-9.
5. Dowton SR and Colten HR (1988) Acute Phase Reactants in Inflammation and Infection. Semin Hematol 25(2): 84-90.
6. Gambino R (1994) C-Reactive Protein (CRP) - How much Proof do we need? Lab Rep 16(11): 83-5
7. Thomas L (2012) Labor und Diagnose – Indikationen und Bewertung von Laborbefunden für die medizinische Diagnostik. TH-Books Verlagsgesellschaft mbH. 8. Auflage. Band 2. S.1279.

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### 16. Referenties

1. Morley JJ, Kushner (1982) Serum C-reactive protein levels in disease. In: Kushner I, Volanakis JE, Gewurz H, eds. C-reactive protein and the plasma protein response to tissue injury. Ann. NY Acad. Sci. 389: 406-417
2. Van Lente F (1982) The Diagnostic Utility of C-Reactive Protein. Hum Path 13(12): 1061-3.
3. Peltola HO (1982) C-reactive protein for rapid monitoring of infections of the central nervous system. Lancet:980-983.
4. Shaw AC (1991) Serum C-Reactive Protein and Neopterin Concentrations in Patients with Viral or Bacterial Infection. J Clin Pathol 44(7): 596-9.

5. Dowton SR and Colten HR (1988) Acute Phase Reactants in Inflammation and Infection. Semin Hematol 25(2): 84-90.
6. Gambino R (1994) C-Reactive Protein (CRP) - How much Proof do we need? Lab Rep 16(11): 83-5
7. Thomas L (2012) Labor und Diagnose – Indikationen und Bewertung von Laborbefunden für die medizinische Diagnostik. TH-Books Verlagsgesellschaft mbH.
8. Auflage. Band 2. S.1279.

Rev.0, 2020-01-28 HB

| Symbol                                                                              | Deutsch                        | English                               | Français                                     | Español                                            | Italiano                                    | Polski                             |
|-------------------------------------------------------------------------------------|--------------------------------|---------------------------------------|----------------------------------------------|----------------------------------------------------|---------------------------------------------|------------------------------------|
|    | CE<br>Konformitätszeichen      | CE marking of<br>conformity           | Conformité aux<br>normes<br>européennes      | Conformidad europea                                | Conformità europea                          | Znak zgodności CE                  |
|    | Gebrauchsanweisung<br>beachten | Consult instructions<br>for use       | Consulter la notice<br>d'utilisation         | Consúltense las<br>instrucciones de uso            | Consultare le<br>istruzioni per l'uso       | Przestrzegać instrukcji<br>obsługi |
|     | In-vitro-Diagnostika           | In-vitro diagnostic<br>medical device | Dispositif médical de<br>diagnostic in vitro | Producto sanitario<br>para diagnóstico<br>in vitro | Dispositivo medico-<br>diagnostico in vitro | Tylko do diagnostyki in<br>vitro   |
|    | Temperaturbegrenzung           | Temperature<br>limitation             | Limites de<br>température                    | Límite de temperatura                              | Limiti di temperatura                       | Temperatura<br>przechowywania      |
|     | Chargenbezeichnung             | Batch code                            | Code du lot                                  | Código de lote                                     | Codice lotto                                | Numer serii                        |
|    | Nicht zur<br>Wiederverwendung  | Do not reuse                          | Ne pas réutiliser                            | No reutilizar                                      | Non riutilizzare                            | Tylko do jednorazowego<br>użytku   |
|    | Verwendbar bis                 | Use by                                | Utiliser jusqu'au                            | Fecha de caducidad                                 | Utilizzare entro                            | Data ważności                      |
|     | Bestellnummer                  | Catalogue Number                      | Référence du<br>catalogue                    | Número de catálogo                                 | Riferimento di<br>Catalogo                  | Numer katalogowy                   |
|   | Hersteller                     | Manufacturer                          | Fabricant                                    | Fabricante                                         | Fabbricante                                 | Producent                          |
|  | Ausreichend für <n><br>Ansätze | Sufficient for<br><n> tests           | Suffisant pour "n"<br>tests                  | Suficiente para <n><br>utilizaciones               | Sufficiente per<br>"n" saggi                | Wystarczający na <n><br>Powtórzeń  |

| Symbol                                                                             | Português                                          | Český                                               | Suomi                                                                | Svenskt                                                     | Nederlands                                        | Dansk                                       | Norsk                                  |
|------------------------------------------------------------------------------------|----------------------------------------------------|-----------------------------------------------------|----------------------------------------------------------------------|-------------------------------------------------------------|---------------------------------------------------|---------------------------------------------|----------------------------------------|
|  | Conformidade<br>com as normas<br>europeias         | CE certifikát                                       | CE-merkitty                                                          | CE-märkning                                                 | CE-markering                                      | CE-mærkning                                 | CE standardisert                       |
|  | Consultar as<br>instruções de<br>utilização        | Viz návod k použití                                 | Katso<br>käyttöohjetta                                               | Läs bruksanvisningen                                        | Raadpleeg de<br>gebruiksaanwijzing                | Se brugsanvisningen                         | Les<br>bruksanvisning<br>nøye          |
|  | Dispositivo médico<br>para diagnóstico in<br>vitro | Diagnostický<br>zdravotnický<br>prostředek in vitro | In vitro -<br>diagnostiikkaan<br>tarkoitettu<br>lääkinnällinen laite | Medicinteknisk<br>produkt avsedd för<br>in vitro-diagnostik | Medisch<br>hulpmiddel voor<br>in-vitrodiagnostiek | Medicinsk udstyr til<br>in vitro-diagnostik | In-vitro diagnostic<br>medisinsk enhet |
|  | Limites de<br>temperatura                          | Teplotní omezení                                    | Lämpötilarajat                                                       | Temperatur-<br>begränsning                                  | Temperatuurlimiet                                 | Temperatur-<br>begrænsning                  | Temperatur<br>begrænsning              |
|  | Código do lote                                     | Kód šarže                                           | Eräkoodi                                                             | Satsnummer                                                  | Code van de partij                                | Batchcode                                   | Merking                                |
|  | Não reutilizar                                     | Pro jednorázové<br>použití                          | Kertakäyttöinen                                                      | Får inte<br>återanvändas                                    | Niet opnieuw<br>gebruiken                         | Må ikke genbruges                           | Må ikke brukes om<br>igjen             |
|  | Prazo de validade                                  | Spotřebujte do                                      | Käytettävä<br>viimeistään                                            | Används före                                                | Houdbaar tot                                      | Udløbsdato                                  | Tidtaking                              |
|  | Número de<br>catálogo                              | Katalogov číslo                                     | Luettelonumero                                                       | Listnummer                                                  | Catalogus nummer                                  | Bestilingsnummer                            | Katalog nummer                         |
|  | Fabricante                                         | Výrobce                                             | Valmistaja                                                           | Tillverkare                                                 | Fabrikant                                         | Fabrikant                                   | Produsent                              |
|  | Suficiente para<br><n> test                        | Dostačuje pro <n><br>testů                          | Lukumäärä <n><br>test                                                | Räcker till <n> test                                        | Voldoende voor<br><n> test                        | Tilstrækkeligt til <n><br>test              | Tilstrækkelig<br>for<n> tester         |

## Our Teams

### Germany:

#### Regensburg

Tel: +49 941 290 10-0  
Fax: +49 941 290 10-50

#### Moers

Tel: +49 2841 99820-0  
Fax: +49 2841 99820-1

#### Austria:

Tel: +49 941 290 10-29  
Free Tel: 0800 291 565  
Fax: +49 290 10-50  
Free Fax: 0800 298 197

#### UK & Ireland:

Tel: +49 941 290 10-18  
Free Tel –UK: 0808 234 1237  
Free Tel – IRE: 1800 555 080  
Fax: +49 290 10-50

#### France:

France Tel: 0800 915 240  
France Fax: 0800 909 493

### Switzerland

Swiss Tel: 0800 564 720  
Swiss Fax: 0800 837 476

### Belgium

Belgium Tel: 0800 718 82  
Belgium Fax: 0800 747 07

### Luxembourg

Lux. Tel: 800 211 16  
Lux. Fax: 800 261 79

### Spain:

Tel: +49 941 290 10-759  
Free Tel: 900 938 315  
Fax: +49 941 290 10-50  
Free Fax: 900 984 992

### Italy:

Tel: +49 941 290 10-34  
Fax: +49 941 290 10-50

### Poland:

Tel: +49 941 290 10-44  
Free Tel: 00 800 491 15 95  
Fax: +49 941 290 10-50  
Free Fax: 00 800 491 15 94

### Portugal:

Tel: +49 941 290 10-735  
Tel. Verde: 800 849 230  
Fax: +49 941 290 10-50  
Fax Verde: 800 849 229

### Netherlands:

Tel: +31 30 75 600  
Free Tel: 0800 0222 890  
Fax: +31 70 30 30 775  
Free Fax: 0800 024 9519

### Nordic countries (Finland, Norway, Sweden, Denmark):

Tel: +31 703075 607  
Free Tel: +45 80 88 87 53  
Tax: +31 703030 775

### Laboratory Diagnostics Team:

Tel: +49 941 290 10-40  
Fax: +49 941 290 10-50

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Carl-Zeiss-Strasse 12 • 47445 Moers • Germany

[www.nal-vonminden.com](http://www.nal-vonminden.com) • [info@nal-vonminden.com](mailto:info@nal-vonminden.com)

Fon: +49 2841 99820-0 • Fax: +49 2841 99820-1