NADAL® FOB75 Test (test cassette)

REF 272942N-20

de	Gebrauchsanweisung	2	cs	Návod k použití	35
en	Instructions for use	7		Symbols	39
fr	Instructions d'utilisation	11		Our Teams	43
es	Instrucciones de uso	16			
it	Istruzioni per l'uso	21			
pl	Sposób użycia	26			
pt	Instruções de Utilização	30			

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15. Referenzen

1. Ferlay, J., et al.; Cancer incidence and mortality worldwide: Sources, methods and major patterns in GLOBOCAN 2012; Int J Cancer (2015); 136(5): E359-86
2. Hardcastle J.D., et al.; Randomised controlled trial of faecal-occult-blood screening for colorectal cancer; Lancet (1996); 348 (9040): 1467-71
3. Kronborg O., et al.; Randomised study of screening for colorectal cancer with faecal-occult-blood test; Lancet (1996); 348 (9040): 1467-71

4. Allison J.E., et al.; Population screening for colorectal cancer means getting FIT: The past, present, and future of colorectal cancer screening using the fecal immunochemical test for hemoglobin (FIT); Gut Liver (2014); 8 (2): 117-30
5. Park D.I.; et al.; Comparison of guaiac-based and quantitative immunochemical fecal occult blood testing in a population at average risk undergoing colorectal cancer screening; Am J Gastroenterol (2010); 105 (9): 2017-25
6. Sultanian R., et al.; The impact of transitioning from guaiac-fecal occult blood testing to fecal immunochemical testing in a Canadian colon cancer screening program (2019) *Journal of the Canadian Association of Gastroenterology*, gwz009, <https://doi.org/10.1093/jcag/gwz009>
7. Konrad G.; Are medication restrictions before FOBT necessary?; Can Fam Physician.(2012); 58(9): 939-948
8. Halloran S.P., et al.; European guidelines for quality assurance in colorectal cancer screening and diagnosis. First Edition DOI <http://dx.doi.org/10.1055/s-0032-1309791> Endoscopy 2012; 44: SE65–SE87
9. Gies A., et. al.; Direct Comparison of Diagnostic Performance of 9 Quantitative Fecal Immunochemical Tests for Colorectal Cancer Screening: Gastroenterology; (2018); 154: 93-104
10. Yamamoto M., Nakama H., Cost-effectiveness analysis of immunochemical occult blood screening for colorectal cancer among three fecal sampling methods: Hepatogastroenterology, 2000 Mar-Apr, 47 (32) 396-399.

Rev. 0, 2022-05-23 OM

1. Intended Use

The NADAL® FOB75 Test is a lateral flow chromatographic immunoassay (immunological faecal occult blood test, iFOBT) for the qualitative detection of human haemoglobin in human faecal specimens. The test is intended for use as an aid in the diagnosis of pathologies of the lower gastrointestinal tract (see section 12 'Limitations'). The test procedure is not automated and requires no special training or qualification. The NADAL® FOB75 Test is designed for professional use only.

2. Introduction and Clinical Significance

Colorectal cancer is one of the most commonly diagnosed types of cancer and one of the leading causes of cancer-related deaths. Screening for occult blood in stool specimens can improve the chances of detecting colorectal cancer at an early stage, thus reducing mortality rates. Previous, commercially available FOB tests were based on the Guaiac method, which requires a special diet prior to testing to avoid inaccurate results caused by, for example, haem or peroxidases found in some foods. For this reason, immunological tests were developed for the specific detection of human haemoglobin in faecal specimens. Several studies have demonstrated that immunochemical methods are superior to the Guaiac method in terms of analytical specificity, diagnostic sensitivity and operability when detecting advanced adenomas and colorectal carcinomas in the lower gastrointestinal tract.

Expected values

Nearly 55% of all colorectal cancer cases occur in more developed regions. Incidence varies worldwide by up to a factor of 10. The highest age-standardised rates can be found in Australia/New Zealand (ASR 44.8 and 32.2 per 100,000 men and women, respectively) and the lowest in West Africa (ASR 4.5 and 3.8 per 100,000 men and women, respectively). Screening for occult ('hidden') blood in stool specimens can increase the chances of the early detection of colorectal cancer and thus reduce mortality rates.

3. Test Principle

The NADAL® FOB75 Test enables the detection of human haemoglobin through the visual interpretation of colour development on the internal test strip. Monoclonal antibodies to human haemoglobin are immobilised in the test line region (T) of the membrane. During the test, the specimen reacts with antibodies to human haemoglobin which are conjugated to coloured particles and precoated onto the conjugate pad of the test cassette. The mixture then migrates along the membrane by capillary action and interacts with the reagents on the membrane. If there is a sufficient amount of haemoglobin in the specimen, a coloured line will develop in the test line region (T) of the membrane. The presence of this coloured line indicates a positive result, whereas its absence indicates a negative result.

The formation of a coloured line in the control line region (C) serves as a procedural control, indicating that the proper volume of specimen has been added and membrane wicking has occurred.

4. Reagents and Materials Supplied

- 20 NADAL® FOB75 test cassettes

- 20 specimen collection tubes containing buffer*
- 1 package insert

*Buffer containing the following preservative: ProClin™ 300: <0.03%.

No hazard labelling for ProClin™ 300 is required according to Regulation (EC) N° 1272/2008 CLP. Concentrations are below the exemption threshold of <0.03%.

5. Additional Materials Required

- Timer
- Optional: Stool collectors
- Optional: Transport bags
- Optional: Patient instructions for faecal specimen collection (Ref. 0705al00524)
- Optional: Specimen collection containers (if an undiluted faecal specimen is required and the dilution step is to be carried out immediately before testing)

6. Storage and Stability

Test kits should be stored at 2-30°C until the indicated expiry date. Test cassettes are stable until the expiry date printed on the foil pouches. Test cassettes must remain in the sealed foil pouches until use. Do not freeze test kits. Do not use tests beyond the expiry date indicated on the packaging. Care should be taken to protect test kit components from contamination. Do not use test kit components if there is evidence of microbial contamination or precipitation. Biological contamination of dispensing equipment, containers or reagents can lead to inaccurate results.

7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Carefully read through the entire instructions for use prior to testing.
- Do not use the test beyond the expiration date indicated on the packaging.
- Do not use test kit components if the primary packaging is damaged.
- Tests are for single use only.
- Do not add specimens to the reaction area (result area).
- In order to avoid contamination, do not touch the reaction area (result area).
- Avoid cross-contamination of specimens by using a new specimen collection tube for each specimen obtained.
- Do not substitute or mix components from different test kits.
- Do not use the buffer if it is discoloured or turbid. Discolouration or turbidity may be a sign of microbial contamination.
- Do not eat, drink or smoke in the area where specimens and test kits are handled.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being assayed.
- Handle all specimens as if they contain infectious agents. Observe established precautions for microbiological risks throughout all procedures and standard guidelines for the appropriate disposal of specimens.
- The test kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals

does not completely guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled in accordance with usual safety precautions (e.g., do not ingest or inhale).

- Temperature can adversely affect test results.
- Used testing materials should be disposed in accordance with local regulations.

8. Specimen Collection and Preparation

The faecal specimen can be transferred into the specimen collection tube containing buffer either directly by the patient or by qualified personnel.

If an undiluted, native faecal specimen is required, the patient can collect the specimen in a dry specimen collection container (not included in the test kit). The faecal specimen is then transferred into the specimen collection tube containing buffer in the doctor's practice or laboratory.

Note: Please instruct the patient on optimal specimen collection, storage and transport.

Notes on specimen collection and storage:

The NADAL® FOB75 Test is intended for use with human faecal specimens only. Patients should not collect specimens during their menstrual period or 3 days before or after it, or if they have bleeding haemorrhoids, blood in their urine or if they experience diarrhoea or strongly forced bowel movements.

Aspirin and other medications taken in excess may cause gastrointestinal irritation, resulting in occult bleeding. If possible, such substances should be discontinued or avoided - following consultation with the doctor - for at least 48 hours prior to testing.

No dietary restrictions are necessary before testing. Alcohol is suspected of irritating the intestines and should be avoided for at least 48 hours prior to testing.

Urine and the excessive dilution of specimens with toilet water can lead to inaccurate test results. It is recommended to collect faecal specimens using a stool collector.


Perform the test as soon as possible after specimen collection. If it is not possible for the patient to bring the specimen collection tube to the doctor immediately, the specimen can be stored in buffer and refrigerated at 2-8°C for a maximum of 3 days. During transportation of the specimen to the doctor, uninterrupted cooling is not required. However, transportation should be carried out quickly, as a slow degradation of haemoglobin may occur at warm temperatures. At room temperature (15-25°C), transportation should not exceed 48 hours. If testing is not carried out in the doctor's practice/laboratory immediately after receiving of the specimen, it can be stored refrigerated there for a maximum of 3 more days. If testing is delayed, faecal specimens in buffer can be frozen at -20°C for up to 2 months. Avoid repeated freeze-thaw cycles.

Haemoglobin breaks down rapidly in undiluted native faecal specimens. Many laboratories therefore recommend that the patient return the collected specimen as soon as possible, preferably on the same day. If a faecal specimen is not to be tested immediately after it is received by the doctor's

practice/laboratory, it should be kept frozen until testing and then tested promptly.

Specimen preparation:


1. Hold the specimen collection tube upright and remove the applicator stick by unscrewing the light-blue cap. Be careful not to spill or splash buffer solution from the collection tube.
2. Collect approx. 50 mg of faeces by inserting the applicator stick into at least 6 different areas of the faecal specimen. Avoid scooping the faecal specimen.
Note: Please be sure to insert the applicator stick 6 times in a row into the faecal specimen. Do NOT put the applicator stick back into the tube in between collecting the 6 specimens. Be careful not to spill fluid from the tube. The quality of the specimen and compliance with these instructions may affect your test result.
3. Put the applicator stick, along with the faecal specimen, back into the specimen collection tube and tightly screw the light-blue cap closed.
4. Shake the specimen collection tube to completely mix the specimen and buffer.

9. Test Procedure

Bring tests, specimens extracted in buffer and/or controls to room temperature (15-30°C) prior to testing.

1. Remove the NADAL® FOB75 test cassette from the sealed foil pouch and use it as soon as possible. The best results will be obtained if the test is performed immediately after opening the foil pouch. Label the test cassette with the patient or control identification.
2. Place the test cassette on a clean and level surface.
3. Shake the specimen collection tube to ensure good specimen dispersion.
4. Unscrew the white cap and break off the tip of the tube using a piece of tissue.
5. Holding the specimen collection tube vertically, transfer 3 drops of the specimen solution into the specimen well (S) of the test cassette.
Avoid trapping air bubbles in the specimen well (S) and do not add any solution to the result area.
6. Start the timer.
7. Wait for the coloured line(s) to appear. Read the test result after exactly 5 minutes. Do not interpret the result after more than 5 minutes.

10. Result Interpretation

Positive:

A coloured line develops in the control line region (C) and another coloured line develops in the test line region (T).

Note: The colour intensity in the test line region (T) may vary depending on the analyte concentration present in the specimen. Any shade of colour in the test line region (T) should be considered positive. Note that this is a qualitative test only and it cannot determine the analyte concentration in the specimen.

Negative:

A coloured line develops in the control line region (C). No line develops in the test line region (T).

Invalid:

The control line (C) fails to appear. Results from any test which has not produced a control line at the specified reading time must be discarded.

Please review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your distributor.

Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for control line failure.

11. Quality Control

The internal procedural control is included in the test cassette: A coloured line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, correct procedural technique and adequate membrane wicking.

Good laboratory practice (GLP) recommends the use of external control materials to ensure proper test kit performance.

12. Limitations

- The NADAL® FOB75 Test is for professional *in-vitro* diagnostic use only. The test should be used for the qualitative detection of haemoglobin in human faecal specimens only.
- Neither the quantitative value nor the rate of increase/decrease in the concentration of human haemoglobin can be determined using this qualitative test.
- Adjustments to the cut-off for FOB tests have been discussed (e.g. according to the patient age and gender) in recent literature (see references 8 and 9). In general, the lower the cut-off, the higher the sensitivity and the lower the specificity. The detection limit for the NADAL® FOB75 Test was chosen so that the test has a high sensitivity. When used for screening purposes, locally recommended guidelines should be followed.
- The accuracy of the test depends on the quality of the specimen. Inaccurate test results may occur due to

improper specimen collection or storage (see section 8 'Specimen Collection and Preparation').

- Other than colorectal bleeding, the presence of blood in the stool may have a number of different causes, such as haemorrhoids, blood in the urine or stomach irritation. Menstrual blood, minor injuries (e.g. from hard stools, rectal medication or digital rectal examinations) or diseases (e.g. diarrhoeal disease, inflammatory bowel disease) may also lead to blood in the stool. Drugs that inhibit blood clotting or irritate the intestinal mucosa are also suspected of exacerbating bleeding tendencies. However, the influence of drugs is controversial, as studies often do not demonstrate any significant differences on the positive predictive value. For a comprehensive review, see reference 7.
- Negative results do not exclude bleeding, as blood may be unevenly distributed in human faecal specimens. In addition, some polyps and colorectal carcinomas can bleed intermittently or not at all. In particular, colorectal polyps in the very early stages often go undetected. Some cancer precursors are detected more reliably using colonoscopy than iFOBT, which is why colonoscopy continues to be regarded as the 'gold standard'. However, as an invasive examination procedure, it poses a greater risk for the patient. Literature demonstrates that multiple measurements over a short period of time significantly increase the diagnostic sensitivity of iFOB tests when compared to a single measurement, with only a slight decrease in diagnostic specificity (See reference 10).
- In rare cases, some untested cross-reacting and interfering substances, as well as faecal components, may interfere with the test and cause inaccurate test results. If clinical signs and symptoms are not consistent with test results, additional testing using other clinical methods (e.g. colonoscopy) is recommended.
- The NADAL® FOB75 Test only detects the presence of human haemoglobin in specimens and should not be used as the sole criterion for a diagnosis.
- The NADAL® FOB75 Test is not suitable for detecting bleeding in the upper gastrointestinal tract, as haemoglobin degrades as it passes through it.
- Not all colorectal bleeding is caused by precancerous or cancerous polyps. As with all diagnostic tests, all results should be interpreted by a physician in conjunction with other available clinical information. Positive test results must be clarified by means of follow-up examinations. Persistent symptoms despite a negative result should be clarified by differential diagnosis.

13. Performance Characteristics

Clinical performance

Diagnostic sensitivity and specificity

The NADAL® FOB75 Test was evaluated using clinical specimens whose status was confirmed using another rapid test.

The results are presented in the following table:

NADAL® FOB75 Test	Another rapid test			
		Positive	Negative	Total
	Positive	80	2	82
	Negative	0	148	148
Total		80	150	230

Diagnostic sensitivity: >99.9% (95.5% - 100%)*

Diagnostic specificity: 98.7% (95.3% - 99.8%)*

Overall agreement: 99.1% (96.9% - 99.9%)*

*95% confidence interval

Analytical performance

Detection limit

The detection limit of the NADAL® FOB75 Test is 75 ng/mL if testing haemoglobin diluted in buffer. In the event of delayed result reading times, specimens containing human haemoglobin at concentrations below 75 ng/mL may also show a positive result.

Measuring range

No adverse effect on T-line formation (prozone effect) was observed when testing specimens containing a haemoglobin concentration as high as 10 µg/mL. At 100 µg/mL a slight decrease in T-line intensity was observed, which is indicative of the beginnings of a prozone effect. However, Hb concentration as high as 1000 µg/mL in buffer were reliably detected as positive. Thus, the measuring range of the test is between 75 ng/mL and 1000 µg/mL (= 3.0 µg/g to 40 mg/g in native faecal specimens if 50 mg faeces are transferred into the buffer). If you suspect a prozone effect, dilute the specimen with the buffer from a new collection tube and repeat the measurement.

Analytical specificity

Cross-reactivity study

The NADAL® FOB75 Test is specific for human haemoglobin and shows no cross-reaction with haemoglobin from bovine, pig, rabbit, chicken, goat and dog blood at concentrations of up to 0.5 mg/mL.

Interference study

Controls without haemoglobin and controls containing a haemoglobin concentration of 75 ng/mL were spiked with the following, potentially interfering substances. The substances showed no interference with the NADAL® FOB75 Test.

Substance	Concentration	Substance	Concentration
Ascorbic acid	0.2 mg/mL	Urea	20 mg/mL
Acetylsalicylic acid	0.2 mg/mL	Glucose	20 mg/mL
Oxalic acid	0.6 mg/mL	Caffeine	0.4 mg/mL
Bilirubin	1 mg/mL	Albumin	20 mg/mL
Uric acid	0.6 mg/mL		

Precision

Repeatability

Repeatability was established by testing 20 replicates of controls at four concentration levels (negative, 40 ng/mL, 75 ng/mL and 1000 ng/mL haemoglobin). Testing was performed by one operator using 3 lots of the NADAL® FOB75

tests at one site on one day. >99% of the specimens were correctly identified (20/20 correct tests per concentration, 95% confidence interval: 95.5% - 100% for each lot). The NADAL® FOB75 Test demonstrated acceptable repeatability.

Reproducibility

Reproducibility was established by testing 5 replicates of controls at four concentration levels (negative, 40 ng/mL, 75 ng/mL and 1000 ng/mL haemoglobin). Testing was performed by 5 operators using 3 independent NADAL® FOB75 test lots at 3 different sites on 5 separate days. >99% of the specimens were correctly identified (375/375 correct tests per concentration, 95% confidence interval: 99.8% - 100%). The NADAL® FOB75 Test demonstrated acceptable reproducibility.

14. Serious incident reporting

In the event of any serious incidents related to the performance of the NADAL® FOB75 Test, please inform nal von minden GmbH and the competent authority immediately. If still possible, do **not** dispose of the used test and the corresponding test kit components.

15. References

1. Ferlay, J., et al.; Cancer incidence and mortality worldwide: Sources, methods and major patterns in GLOBOCAN 2012; Int J Cancer (2015); 136(5): E359-86
2. Hardcastle J.D., et al.; Randomised controlled trial of faecal-occult-blood screening for colorectal cancer; Lancet (1996); 348 (9040): 1467-71
3. Kronborg O., et al.; Randomised study of screening for colorectal cancer with faecal-occult-blood test; Lancet (1996); 348 (9040): 1467-71
4. Allison J.E., et al.; Population screening for colorectal cancer means getting FIT: The past, present, and future of colorectal cancer screening using the fecal immunochemical test for hemoglobin (FIT); Gut Liver (2014); 8 (2): 117-30
5. Park D.I., et al.; Comparison of guaiac-based and quantitative immunochemical fecal occult blood testing in a population at average risk undergoing colorectal cancer screening; Am J Gastroenterol (2010); 105 (9): 2017-25
6. Sultanian R., et al.; The impact of transitioning from guaiac-fecal occult blood testing to fecal immunochemical testing in a Canadian colon cancer screening program (2019) *Journal of the Canadian Association of Gastroenterology*, gwz009, <https://doi.org/10.1093/jcag/gwz009>
7. Konrad G.; Are medication restrictions before FOBT necessary?; Can Fam Physician.(2012); 58(9): 939-948
8. Halloran S.P., et al.; European guidelines for quality assurance in colorectal cancer screening and diagnosis. First Edition DOI <http://dx.doi.org/10.1055/s-0032-1309791> Endoscopy 2012; 44: SE65–SE87
9. Gies A., et. al.; Direct Comparison of Diagnostic Performance of 9 Quantitative Fecal Immunochemical Tests for Colorectal Cancer Screening; Gastroenterology; (2018); 154: 93-104
10. Yamamoto M., Nakama H.; Cost-effectiveness analysis of immunochemical occult blood screening for colorectal cancer among three fecal sampling methods: Hepatogastroenterology, 2000 Mar-Apr, 47 (32) 396-399.

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15. Bibliographie

- Ferlay, J., et al.; Cancer incidence and mortality worldwide: Sources, methods and major patterns in GLOBOCAN 2012; Int J Cancer (2015); 136(5): E359-86
- Hardcastle J.D., et al.; Randomised controlled trial of faecal-occult-blood screening for colorectal cancer; Lancet (1996); 348 (9040): 1467-71
- Kronborg O., et al.; Randomised study of screening for colorectal cancer with faecal-occult-blood test; Lancet (1996); 348 (9040): 1467-71
- Allison J.E., et al.; Population screening for colorectal cancer means getting FIT: The past, present, and future of colorectal cancer screening using the fecal immunochemical test for hemoglobin (FIT); Gut Liver (2014); 8 (2): 117-30

5. Park D.I., et al.; Comparison of guaiac-based and quantitative immunochemical fecal occult blood testing in a population at average risk undergoing colorectal cancer screening; Am J Gastroenterol (2010); 105 (9): 2017-25
6. Sultanian R., et al.; The impact of transitioning from guaiac-fecal occult blood testing to fecal immunochemical testing in a Canadian colon cancer screening program (2019) *Journal of the Canadian Association of Gastroenterology*, gwz009, <https://doi.org/10.1093/jcag/gwz009>
7. Konrad G.; Are medication restrictions before FOBT necessary?; Can Fam Physician.(2012); 58(9): 939-948
8. Halloran S.P., et al.; European guidelines for quality assurance in colorectal cancer screening and diagnosis. First Edition DOI <http://dx.doi.org/10.1055/s-0032-1309791> Endoscopy 2012; 44: SE65–SE87
9. Gies A., et. al.; Direct Comparison of Diagnostic Performance of 9 Quantitative Fecal Immunochemical Tests for Colorectal Cancer Screening: Gastroenterology; (2018); 154: 93-104
10. Yamamoto M., Nakama H., Cost-effectiveness analysis of immunochemical occult blood screening for colorectal cancer among three fecal sampling methods: Hepatogastroenterology, 2000 Mar-Apr, 47 (32) 396-399.

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15. Referencias

1. Ferlay, J., et al.; Cancer incidence and mortality worldwide: Sources, methods and major patterns in GLOBOCAN 2012; Int J Cancer (2015); 136(5): E359-86
2. Hardcastle J.D., et al.; Randomised controlled trial of faecal-occult-blood screening for colorectal cancer; Lancet (1996); 348 (9040): 1467-71
3. Kronborg O., et al.; Randomised study of screening for colorectal cancer with faecal-occult-blood test; Lancet (1996); 348 (9040): 1467-71
4. Allison J.E., et al.; Population screening for colorectal cancer means getting FIT: The past, present, and future of colorectal cancer screening using the fecal immunochemical test for hemoglobin (FIT); Gut Liver (2014); 8 (2): 117-30
5. Park D.I., et al.; Comparison of guaiac-based and quantitative immunochemical fecal occult blood testing in a population at average risk undergoing colorectal cancer screening; Am J Gastroenterol (2010); 105 (9): 2017-25
6. Sultanian R., et al.; The impact of transitioning from guaiac-fecal occult blood testing to fecal immunochemical testing in a Canadian colon cancer screening program (2019) Journal of the Canadian Association of Gastroenterology, gwz009, <https://doi.org/10.1093/jcag/gwz009>
7. Konrad G.; Are medication restrictions before FOBT necessary?; Can Fam Physician.(2012); 58(9): 939-948

8. Halloran S.P., et al.; European guidelines for quality assurance in colorectal cancer screening and diagnosis. First Edition DOI <http://dx.doi.org/10.1055/s-0032-1309791> Endoscopy 2012; 44: SE65–SE87
9. Gies A., et. al.; Direct Comparison of Diagnostic Performance of 9 Quantitative Fecal Immunochemical Tests for Colorectal Cancer Screening: Gastroenterology; (2018); 154: 93-104
10. Yamamoto M., Nakama H., Cost-effectiveness analysis of immunochemical occult blood screening for colorectal cancer among three fecal sampling methods: Hepatogastroenterology, 2000 Mar-Apr, 47 (32) 396-399.

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15. Bibliografia

1. Ferlay, J., et al.; Cancer incidence and mortality worldwide: Sources, methods and major patterns in GLOBOCAN 2012; Int J Cancer (2015); 136(5): E359-86
2. Hardcastle J.D., et al.; Randomised controlled trial of faecal-occult-blood screening for colorectal cancer; Lancet (1996); 348 (9040): 1467-71
3. Kronborg O., et al.; Randomised study of screening for colorectal cancer with faecal-occult-blood test; Lancet (1996); 348 (9040): 1467-71
4. Allison J.E., et al.; Population screening for colorectal cancer means getting FIT: The past, present, and future of colorectal cancer screening using the fecal immunochemical test for hemoglobin (FIT); Gut Liver (2014); 8 (2): 117-30

5. Park D.I., et al.; Comparison of guaiac-based and quantitative immunochemical fecal occult blood testing in a population at average risk undergoing colorectal cancer screening; *Am J Gastroenterol* (2010); 105 (9): 2017-25
6. Sultanian R., et al.; The impact of transitioning from guaiac-fecal occult blood testing to fecal immunochemical testing in a Canadian colon cancer screening program (2019) *Journal of the Canadian Association of Gastroenterology*, gwz009, <https://doi.org/10.1093/jcag/gwz009>
7. Konrad G.; Are medication restrictions before FOBT necessary?; *Can Fam Physician*.(2012); 58(9): 939-948
8. Halloran S.P., et al.; European guidelines for quality assurance in colorectal cancer screening and diagnosis. First Edition DOI <http://dx.doi.org/10.1055/s-0032-1309791> *Endoscopy* 2012; 44: SE65–SE87
9. Gies A., et. al.; Direct Comparison of Diagnostic Performance of 9 Quantitative Fecal Immunochemical Tests for Colorectal Cancer Screening: *Gastroenterology*; (2018); 154: 93-104
10. Yamamoto M., Nakama H., Cost-effectiveness analysis of immunochemical occult blood screening for colorectal cancer among three fecal sampling methods: *Hepatogastroenterology*, 2000 Mar-Apr, 47 (32) 396-399.

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15. Bibliografia

1. Ferlay, J., et al.; Cancer incidence and mortality worldwide: Sources, methods and major patterns in GLOBOCAN 2012; Int J Cancer (2015); 136(5): E359-86
2. Hardcastle J.D., et al.; Randomised controlled trial of faecal-occult-blood screening for colorectal cancer; Lancet (1996); 348 (9040): 1467-71
3. Kronborg O., et al.; Randomised study of screening for colorectal cancer with faecal-occult-blood test; Lancet (1996); 348 (9040): 1467-71
4. Allison J.E., et al.; Population screening for colorectal cancer means getting FIT: The past, present, and future of colorectal cancer screening using the fecal immunochemical test for hemoglobin (FIT); Gut Liver (2014); 8 (2): 117-30
5. Park D.I., et al.; Comparison of guaiac-based and quantitative immunochemical fecal occult blood testing in a population at average risk undergoing colorectal cancer screening; Am J Gastroenterol (2010); 105 (9): 2017-25
6. Sultanian R., et al.; The impact of transitioning from guaiac-fecal occult blood testing to fecal immunochemical testing in a Canadian colon cancer screening program (2019) Journal of the Canadian Association of Gastroenterology, gwz009, <https://doi.org/10.1093/jcag/gwz009>
7. Konrad G.; Are medication restrictions before FOBT necessary?; Can Fam Physician (2012); 58(9): 939-948
8. Halloran S.P., et al.; European guidelines for quality assurance in colorectal cancer screening and diagnosis. First Edition DOI <http://dx.doi.org/10.1055/s-0032-1309791> Endoscopy 2012; 44: SE65-SE87
9. Gies A., et al.; Direct Comparison of Diagnostic Performance of 9 Quantitative Fecal Immunochemical Tests for Colorectal Cancer Screening; Gastroenterology; (2018); 154: 93-104
10. Yamamoto M., Nakama H.; Cost-effectiveness analysis of immunochemical occult blood screening for colorectal cancer among three fecal sampling methods: Hepatogastroenterology, 2000 Mar-Apr, 47 (32) 396-399.

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15. Reference

1. Ferlay, J., et al.; Cancer incidence and mortality worldwide: Sources, methods and major patterns in GLOBOCAN 2012; Int J Cancer (2015); 136(5): E359-86
2. Hardcastle J.D., et al.; Randomised controlled trial of faecal-occult-blood screening for colorectal cancer; Lancet (1996); 348 (9040): 1467-71
3. Kronborg O., et al.; Randomised study of screening for colorectal cancer with faecal-occult-blood test; Lancet (1996); 348 (9040): 1467-71
4. Allison J.E., et al.; Population screening for colorectal cancer means getting FIT: The past, present, and future of colorectal cancer screening using the fecal immunochemical test for hemoglobin (FIT); Gut Liver (2014); 8 (2): 117-30
5. Park D.I., et al.; Comparison of guaiac-based and quantitative immunochemical fecal occult blood testing in a population at average risk undergoing colorectal cancer screening; Am J Gastroenterol (2010); 105 (9): 2017-25
6. Sultanian R., et al.; The impact of transitioning from guaiac-fecal occult blood testing to fecal immunochemical testing in a Canadian colon cancer screening program (2019) *Journal of the Canadian Association of Gastroenterology*, gwz009, <https://doi.org/10.1093/jcag/gwz009>
7. Konrad G.; Are medication restrictions before FOBT necessary?; Can Fam Physician (2012); 58(9): 939-948
8. Halloran S.P., et al.; European guidelines for quality assurance in colorectal cancer screening and diagnosis. First Edition DOI <http://dx.doi.org/10.1055/s-0032-1309791> Endoscopy 2012; 44: SE65-SE87
9. Gies A., et al.; Direct Comparison of Diagnostic Performance of 9 Quantitative Fecal Immunochemical Tests for Colorectal Cancer Screening; Gastroenterology; (2018); 154: 93-104
10. Yamamoto M., Nakama H.; Cost-effectiveness analysis of immunochemical occult blood screening for colorectal cancer among three fecal sampling methods: Hepatogastroenterology, 2000 Mar-Apr, 47 (32) 396-399.

Rev. 0, 2022-05-23 JV

Symbol	Deutsch	English	Français	Español	Italiano	Polski
	CE Konformitätszeichen	CE marking of conformity	Conformité aux normes européennes	Conformidad europea	Conformità europea	Znak zgodności CE
	Gebrauchsanweisung beachten	Consult instructions for use	Consulter la notice d'utilisation	Consúltense las instrucciones de uso	Consultare le istruzioni per l'uso	Przestrzegać instrukcji obsługi
	<i>in-vitro</i> -Diagnostika	<i>in-vitro</i> diagnostic medical device	Dispositif médical de diagnostic <i>in vitro</i>	Producto sanitario para diagnóstico <i>in-vitro</i>	Dispositivo medico- diagnostico <i>in-vitro</i>	Tylko do diagnostyki <i>in-vitro</i>
	Temperaturbegrenzung	Temperature limitation	Limites de température	Límite de temperatura	Limiti di temperatura	Temperatura przechowywania
	Chargenbezeichnung	Batch code	Code du lot	Código de lote	Codice lotto	Numer serii
	Nicht zur Wiederverwendung	Do not reuse	Ne pas réutiliser	No reutilizar	Non riutilizzare	Tylko do jednorazowego użytku
	Verwendbar bis	Use by	Utiliser jusqu'au	Fecha de caducidad	Utilizzare entro	Data ważności
	Bestellnummer	Catalogue Number	Référence du catalogue	Número de catálogo	Riferimento di Catalogo	Numer katalogowy
	Hersteller	Manufacturer	Fabricant	Fabricante	Fabbricante	Producent
	Ausreichend für <n> Ansätze	Sufficient for <n> tests	Suffisant pour "n" tests	Suficiente para <n> utilizaciones	Sufficiente per "n" saggi	Wystarczający na <n> Powtórzeń

Symbol	Português	Český	Suomi	Svenskt	Nederlands	Dansk	Norsk
	Conformidade com as normas europeias	CE certifikát	CE-merkitty	CE-märkning	CE-markering	CE-mærkning	CE standardisert
	Consultar as instruções de utilização	Viz návod k použití	Katso käyttöohjetta	Läs bruksanvisningen	Raadpleeg de gebruiksaanwijzing	Se brugsanvisningen	Les bruksanvisning nøye
	Dispositivo médico para diagnóstico <i>in-vitro</i>	Diagnostický zdravotnický prostředek <i>in-vitro</i>	<i>in-vitro</i> - diagnostiikkaan tarkoitettu lääkinnällinen laite	Medicinteknisk produkt avsedd för <i>in-vitro</i> -diagnostik	Medisch hulpmiddel voor <i>in-vitro</i> diagnostiek	Medicinsk udstyr til <i>in-vitro</i> -diagnostik	<i>in-vitro</i> diagnostic medisinsk enhet
	Limites de temperatura	Teplotní omezení	Lämpötilarajat	Temperatur- begränsning	Temperatuurlimiet	Temperatur- begrænsning	Temperatur begrænsning
	Código do lote	Kód šarže	Eräkoodi	Satsnummer	Code van de partij	Batchkode	Merking
	Não reutilizar	Pro jednorázové použití	Kertakäyttöinen	Får inte återanvändas	Niet opnieuw gebruiken	Må ikke genbruges	Må ikke brukes om igjen
	Prazo de validade	Spotřebujte do	Käytettävä viimeistään	Används före	Houdbaar tot	Udløbsdato	Tidtaking
	Número de catálogo	Katalogov číslo	Luettelonumero	Listnummer	Catalogus nummer	Bestilingsnummer	Katalog nummer
	Fabricante	Výrobce	Valmistaja	Tillverkare	Fabrikant	Fabrikant	Produsent
	Suficiente para <n> test	Dostačuje pro <n> testů	Lukumäärä <n> test	Räcker till <n> test	Voldoende voor <n> test	Tilstrækkeligt til <n> test	Tilstrækkelig for<n> tester

Our Teams**Germany:****Regensburg**

Tel: +49 941 290 10-0
Fax: +49 941 290 10-50

Moers

Tel: +49 2841 99820-0
Fax: +49 2841 99820-1

Austria:

Tel: +49 941 290 10-29
Free Tel: 0800 291 565
Fax: +49 290 10-50
Free Fax: 0800 298 197

UK & Ireland:

Tel: +49 941 290 10-18
Free Tel –UK: 0808 234 1237
Free Tel – IRE: 1800 555 080
Fax: +49 290 10-50

France:

France Tel: 0800 915 240
France Fax: 0800 909 493

Switzerland

Swiss Tel: 0800 564 720
Swiss Fax: 0800 837 476

Belgium

Belgium Tel: 0800 718 82
Belgium Fax: 0800 747 07

Luxembourg

Lux, Tel: 800 211 16
Lux, Fax: 800 261 79

Spain:

Tel: +49 941 290 10-759
Free Tel: 900 938 315
Fax: +49 941 290 10-50
Free Fax: 900 984 992

Italy:

Tel: +49 941 290 10-34
Fax: +49 941 290 10-50

Poland:

Tel: +49 941 290 10-44
Free Tel: 00 800 491 15 95
Fax: +49 941 290 10-50
Free Fax: 00 800 491 15 94

Portugal:

Tel: +49 941 290 10-735
Tel, Verde: 800 849 230
Fax: +49 941 290 10-50
Fax Verde: 800 849 229

Netherlands:

Tel: +31 30 75 600
Free Tel: 0800 0222 890
Fax: +31 70 30 30 775
Free Fax: 0800 024 9519

Nordic countries:**Denmark**

Tel: +31 703075 605
Free Tel: 808 887 53

Finland

Tel: +31 703075 606
Free Tel: 0800 918 263
Free Fax: 0800 918 262

Norway

Tel: +31 703075 605
Free Tel: 800 16 731

Sweden

Tel: +31 703075 605
Free Tel: 020 79 09 06

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Carl-Zeiss-Strasse 12 • 47445 Moers • Germany

www.nal-vonminden.com • info@nal-vonminden.com

Tel: +49 2841 99820-0 • Fax: +49 2841 99820-1